









Product data sheet – ReadyCell Procedure Pack

The ReadyCell Procedure Pack (RCP-100) consists of a package of accessories to be associated with the Sepax® ReadyCell protocol for the preparation of cellular therapy products. In combination with the Sepax® 2 processing platform, it offers an all-encompassing solution from bone marrow harvesting to safe and fast concentration of stem cells and reinjection back into the patient.

Components of RCP-100 in order of procedure steps:		
	2 x aspiration needles	Needles for bone marrow harvesting; 11 gauge x 100 mm, with ergonomic handle, luer lock connection for syringe and remover guide with indicators.
	3 x 30 ml luer-lock syringes	Syringes to be used for bone marrow aspiration; include a patented system to block the piston for a better secured aspiration.
	4 x 10 ml luer-lock syringes	Alternative syringes for bone marrow aspiration in case a smaller diameter of syringe is needed.
	Collection bag RCA-100 with integrated filter	Collection bag for bone marrow harvest; includes filter for addition of anti-coagulant and sampling port. To be connected to CS-470.1 single-use kit with included spike or sterile connection.
	2 x 5 ml luer-lock syringes	Syringes for initial sampling of bone marrow.
	CS-470.1 single-use kit	Sepax single-use processing kit specifically designed for the ReadyCell protocol.

	2 x cannulas	Cannulas for stem cell retrieval from CS-470.1 vial; 14 gauge x 80 mm
	20 ml luer-lock syringe	Syringe for stem cell retrieval from CS-470.1 vial.

Storage requirements:

The ReadyCell Procedure Pack must be stored in a clean and dry place, free from biological contaminants and chemical vapors. It must only be used/ stored under the following environmental conditions:

Temperature	Operation	During Transport
Relative humidity	+7°C to +27°C 30% to 75%, non-condensing	4°C to 40°C 20% to 75%, non-condensing
<ul style="list-style-type: none"> • Shelf life • Sterility • Material 	18 months from manufacturing All components of the RCP-100 are provided sterile by their supplier All material are latex free	

Order information:

RCP-100	Product #10035
Packaging	Individual alveolar box (box dimensions: 37x44.3x7.5 cm/ 14.5x17.4x2.9 in)
Minimum order quantity	Box of 3 kits (box dimension: 38.9x50x26.1 cm/ 15.3x19.7x10.3 in)

All products included in the RCP-100 are CE marked.

Head Office
Biosafe Group SA

Eysins, Switzerland
T: +41 22 365 27 27
E: info@biosafe.ch

North America Office
Biosafe America, Inc.

Houston, USA
T: +1 713 936 2615
E: info@biosafeamerica.com

Latin America office
Biosafe Latin America

São Paulo, Brazil
T : +55 11 2655 1731
E: info@biosafelatinamerica.com

Asia office
Biosafe Asia-Pacific Ltd.

Hong Kong
T: +852 2956 7500
E : info@biosafeasiapacific.com

China office
Biosafe Medical Device Int. Trading (Shanghai) Co. Ltd

Shanghai, China
T: +86 21 6137 3214
E : info@biosafechina.com